


# Breaking down barriers + Progressing wounds


IODOSORB<sup>®</sup> is indicated for use in chronic and infected wounds where bioburden is a barrier to healing.

Prepare the wound to progress to closure.




**IODOSORB** <sup>◇</sup> unique antimicrobial dressing made of cadexomer micro-beads: spherical starch structures loaded with 0.9% elemental iodine.

IODOSORB changes from brown to yellow/white as the iodine is exhausted<sup>26-30</sup>



### Ordering codes for IODOSORB

IODOSORB Ointment 	Size	Product #	Pcs/Pkg
	10g	66060630	4
	20g	66060631	2
IODOSORB Dressing 	40g	66060632	1
	5g (4cm x 6cm)	66060633	5
	10g (6cm x 8cm)	66060634	3
IODOSORB Powder 	17g (8cm x 10cm)	66060635	2
	3g (7 x 3G)	66001286	7

Advanced Wound Management  
Smith & Nephew Inc.  
2280 Argentia Road  
Mississauga, ON L5N 6H8,  
Canada

[www.smith-nephew.com](http://www.smith-nephew.com)

T1-800-463-7439  
F1-800-671-9140

\*Trademark of Smith & Nephew  
All Trademarks acknowledged  
© 2019 Smith & Nephew  
CA32707 04/20

References: 1. Skog E, et al. A randomized trial comparing cadexomer iodine and standard treatment in the out-patient management of chronic venous ulcers. *Br. J. Dermatol.* 109, 77–83 (1983). 2. Trøeng T, Skog E, Arnesjö B, Gjores JE, Bergljung L, Gundersen JEA. A randomised multicentre trial to compare the efficacy of cadexomer iodine and standard treatment in the management of chronic venous ulcers in out patients. In *Cadexomer Iodine*, (eds. Fox, J. & Fisher, H.) 43–50 (Schattauer Verlag, 1983). 3. Ormiston MC, Fox JA. Controlled trial of Iodosorb. *Br. Med. J. (Clin. Res. Ed)*. 291, 308–310 (1985) 4. Hansson C, et al. The effects of cadexomer iodine paste in the treatment of venous leg ulcers compared with hydrocolloid dressing and paraffin gauze dressing. *Int. J. Dermatol.* 37, 390–396 (1998). 5. Lindsay G, Latta D, Lyons KGB. A study in general practice of the efficacy of cadexomer iodine in venous leg ulcers treated on alternate days. *Acta Ther.* 12, 141–148 (1986). 6. Johnson A. A combative healer with no ill-effect. Iodosorb in the treatment of infected wounds. *Prof. Nurse* 7, 60, 62, 64 (1991). 7. Harcup JW, Saul PA. A study of the effect of cadexomer iodine in the treatment of venous leg ulcers. *Br. J. Clin. Pract.* 40, 360–4 (1986). 8. Hesler B. IODOSORB™ Powder wound model testing. 3 day cavity testing with ALLEVYN™ GENTLE Border secondary dressing using horse serum at a flow rate modelling that of a moderately exuding wound. Smith&Nephew Data on File # DS.17.249.R. (2017). 9. Hesler B. IODOSORB Ointment/Gel wound model testing. 3 day cavity testing with ALLEVYN gentle border secondary dressing using horse serum at a flow rate modelling that of a moderately exuding wound. Smith&Nephew Data on File DS.17.572.R. (2017). 10. Hesler B. IODOSORB IODOFLEX™ dressing wound model testing. 3 day cavity testing with ALLEVYN gentle border secondary dressing using horse serum at a flow rate modelling that of a moderately exuding wound. Smith&Nephew Data on File # DS.17.247.R. (2017). 11. Fitzgerald, D. J. et al. *Wound Repair Regen.* 25, 13–24 (2017). 12. Akiyama, H., Oono, T., Saito, M. & Iwatsuki, K. *J. Dermatol.* 31, 529–34 (2004). 13. Roche, E. et al. *Int. Wound J.* In press, 1–10 (2019). 14. Oates, J. L., Phillips, C. D., Wolcott, R. & Woodmansey, E. in *SAWC Atlanta 1* (2016). 15. Phillips, P. L. et al. *Int. Wound J.* 469–483 (2013). 16. Hill, K. E. et al. *J. Antimicrob. Chemother.* 65, 1195–1206 (2010). 17. Kiamco, M. M., Atci, E., Mohamed, A., Call, D. R. & Beyenal, H. *Appl. Environ. Microbiol.* 1–33 (2017) 18. Forest E. Antimicrobial activity of IODOSORB range against a broad spectrum of wound pathogens. Smith&Nephew Data on File #1801001. (2018). 19. Harrow J. A comparison of the antimicrobial activity of a cadexomer iodine dressing and a povidone iodine dressing. 6763/IODOSORB/TECHMON/GLOBAL/0404. (2009). 20. Sundberg J, Pharm MS, Meller R, Se B, Pharm MS. A retrospective Review of the Use of Cadexomer Iodine in the Treatment of Chronic Wounds. *Wounds A Compend. Clin. Res. Pract.* 9, 68–86 (1997). 21. Akiyama H, Oono T, Saito M, Iwatsuki K. Assessment of cadexomer iodine against *Staphylococcus aureus* biofilm in vivo and in vitro using confocal laser scanning microscopy. *J. Dermatol.* 31, 529–34 (2004). 22. Zhou LH, Nahm WK, Badiavas E, Yufit T, Falanga V. Slow release iodine preparation and wound healing: in vitro effects consistent with lack of in vivo toxicity in human chronic wounds. *Br. J. Dermatol.* 146, 365–74 (2002). 23. Holloway GA, Johansen KH, Barnes RW, Pierce GE. Multicenter trial of cadexomer iodine to treat venous stasis ulcer. *West. J. Med.* 151, 35–38 (1989). 24. Moberg S, Hoffman MOL, Holst A. A Randomized Trial of Cadexomer in Decubitus Ulcers Iodine. 25. Laudanska H, Gustavson B. In-patient treatment of chronic varicose venous ulcers. A randomized trial of cadexomer iodine versus standard dressings. *J. Int. Med.* Res. 16, 428–35 (1988) 26. Hesler B. Release of iodine from IODOSORB/IODOFLEX dressing into 0.9% saline solution over a three day period. Smith&Nephew Data on File # DS.18.026.R. (2018). 27. Hesler B. Release of iodine from IODOSORB powder into 0.9% saline solution over a three day period. Smith&Nephew Data on File # DS.18.024.R. (2018). 28. Hesler B. Release of iodine from IODOSORB ointment into 0.9% saline solution over a three day period. Smith&Nephew Data on File # DS.18.025.R. (2018). 29. Moberg S, Hoffman L, Grennert ML, Holst A. A randomized trial of cadexomer iodine in decubitus ulcers. *J. Am. Geriatr. Soc.* 31, 462–465 (1983). 30. Hesler B. IODOSORB Powder wound model testing to investigate colour loss of the product over the course of a 3 day period. Data on file # DS.17.594.R. (2017). 32. Hesler B. (2017) Method transfer report Smith+Nephew Data on File # DS.17.210.R.